## ASP™

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## US

### EN: Indication for use

Medium for oocyte retrieval and rinsing (follicle flushing).

### Product Description

ASP™ is a bicarbonate and HEPES buffered medium containing heparin and gentamicin as an antibacterial agent. For use after equilibration at +37° and 5% CO2 atmosphere.

### Storage instructions and stability

ASP™ is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

### Directions for use

Equilibrate ASP™ at +37°С and 5% CO2. Rinse the aspiration needle lumen and tubing using ASP™ medium and discard the rinsed fluid.

The follicles may be aspirated individually or collectively. The aspirated follicles should be collected in a sterile, non-toxic tissue culture tube. It is recommended to keep ASP™ at +37°С during the procedure by using thermostat-warming blocks, and in a sealed container to retain pH. The follicle aspirates should be collected in a sterile, tissue culture appropriate vessel and should be examined by the laboratory immediately. If they cannot be examined straight away, the vessels should be tightly sealed and kept at +37°С.

### Specifications

<p>Mouse Embryo Assay (1-cell) [% expanded blastocyst (within 96 hours)]</p> <p>Bacterial endotoxins (LAL assay) [EU/mL]</p> <p>LOT specific test results are available on the Certificate of Analysis provided with each delivery.</p>	<p>SAL 10<sup>+</sup></p> <p>≥ 80</p> <p>&lt; 0.25</p>
---	--

### Precautions

Discard product if bottle integrity is compromised. Do not use ASP™ if it appears cloudy.

ASP™ contains heparin and gentamicin. Do not use in patients with known hypersensitivity/allergy to any of the components.

Re-use may result in microbiological contamination and/or property changes in the product.

To avoid contamination VitroLife strongly recommends that media should be opened and used only with aseptic technique.

The risks of reproductive toxicity and developmental toxicity for IVF media, including VitroLife’s IVF media, have not been determined and are under investigation.

Not to be used as an injectable product other than follicle flushing during oocyte retrieval.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

### Description of ISO Symbols

	<p>Sterilized using aseptic processing techniques</p>
	<p>Temperature limit</p>
	<p>Do not re-use, discard after use</p>
	<p>Date of manufacture</p>
	<p>Use-by date</p>
	<p>Caution: Consult accompanying documents</p>
	<p>Catalog number</p>
	<p>Batch code</p>
	<p>CE Mark (Conformité Européen)</p>

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