# HSA-Solution™

## US

### EN: Indication for use

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## EUROPE

### EN: Indication for use

HSA-solution™ contains human serum albumin solution (100 mg/mL) and is intended for use in assisted reproductive procedures which include gamete and embryo manipulation. These procedures include the use of HSA-solution™ as a supplement for culture medium.

### Product Description

HSA-solution™ is antibiotic-free.

### Storage instructions and stability

Store dark at +2 to +8 °C.

HSA-solution™ is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis. Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

### Directions for use

G-Series PLUS media, G-GAMET™ and EmbryoGlue® do not require albumin supplementation. All other G-Series culture media, not designated PLUS, come protein-free and need to be supplemented either with G-MM™ or HSA-solution™ at the appropriate concentration as listed in the table below. When selecting protein supplementation it should be noted that human serum albumin (HSA) is a blood derived substance and may contain human pathogenic agents including those not yet known or identified. Thus the risk of transmission of such infectious agents cannot be completely eliminated when using HSA. G-Series media with the exception of G-IVF™, G-FreezeKItBlast™ and G-ThawKIt Blast™ shall have 5 % of either G-MM™ or HSA-solution™ added (0.5 mL added to 9.5 mL). The final concentration of G-MM™ will be 2.5 mg/mL and the final concentration of HSA will be 5 mg/mL.

G-IVF™ shall have 10 % of either G-MM™ or HSA-solution™ added (1.0 mL added to 9.0 mL). The final concentration of G-MM™ will be 5 mg/mL and the final concentration of HSA will be 10 mg/mL. Media to be supplemented should be aliquoted into sterile non-toxic tissue culture grade tubes or flasks. All containers should be pre-rinsed using G-RINSE™ in order to ensure there is no particulate matter or toxic substances residing.

### Supplementation of G-MOPS™, G-1™, G-2™, G-PGD™:

| Medium | HSA-solution™ | Final       |
|--------|---------------|-------------|
| [mL]   | volume [mL]   | volume [mL] |
| 9.0    | 0.5           | 10.0        |
| 19.0   | 1.0           | 20.0        |
| 28.5   | 1.5           | 30.0        |
| 38.0   | 2.0           | 40.0        |
| 47.5   | 2.5           | 50.0        |
| 54.0   | 3.0           | 60.0        |
| 63.0   | 3.5           | 70.0        |
| 72.0   | 4.0           | 80.0        |
| 81.0   | 4.5           | 90.0        |
| 95.5   | 5.0           | 100.0       |

The cryopreservation solutions can be supplemented directly into the dishes. Simply place 900 µl of each cryopreservation solution into separate wells + 100 µl of G-MM™ or 100 µl HSA-solution™. It is recommended only to supplement the volume of medium expected to be used in one day.

### Specifications

|                                          |                     |
|------------------------------------------|---------------------|
| Sterile filtered                         | SAL 10 <sup>3</sup> |
| Mouse Embryo Assay (1-cell)              |                     |
| % expanded blastocysts (within 96 hours) | ≥ 80                |
| Bacterial endotoxins (LAL assay) [EU/mL] | < 0.5               |

LOT specific test results are available on the Certificate of Analysis provided with each delivery.

### Precautions

Discard product if bottle integrity is compromised. Do not use HSA-solution™ if it appears cloudy.

HSA-solution™ contains human serum albumin.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested for antibodies to HIV, HBc, HCV, and HTLV III and when not-reactive for HbsAg, HCV RNA and HIV-1 RNA and syphilis. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

Re-use may result in microbiological contamination and/or property changes in the product.

To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique.

The risk of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife's IVF media, have not been determined and are uncertain.

Not for injection.

Caution: Federal (US) law restricts this device to sale by or on behalf of a physician.

(100 mg/mL) and is intended for use in assisted reproductive procedures which include gamete and embryo manipulation. These procedures include the use of HSA-solution™ as a supplement for culture medium.

### Product Description

HSA-solution™ is antibiotic-free.

### Storage instructions and stability

Store dark at +2 to +8 °C.

HSA-solution™ is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis. Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

### Directions for use

G-Series PLUS media, G-GAMET™ and EmbryoGlue® do not require albumin supplementation. All other G-Series culture media, not designated PLUS, come protein-free and need to be supplemented either with G-MM™ or HSA-solution™ at the appropriate concentration as listed in the table below. When selecting protein supplementation it should be noted that human serum albumin (HSA) is a blood derived substance and may contain human pathogenic agents including those not yet known or identified. Thus the risk of transmission of such infectious agents cannot be completely eliminated when using HSA. G-Series media with the exception of G-IV™, G-FreezeKItBlast™ and G-ThawKIt Blast™ shall have 5 % of either G-MM™ or HSA-solution™ added (0.5 mL added to 9.5 mL). The final concentration of G-MM™ will be 2.5 mg/mL and the final concentration of HSA will be 5 mg/mL.

G-IVF™ shall have 10 % of either G-MM™ or HSA-solution™ added (1.0 mL added to 9.0 mL). The final concentration of G-MM™ will be 5 mg/mL and the final concentration of HSA will be 10 mg/mL. Media to be supplemented should be aliquoted into sterile non-toxic tissue culture grade tubes or flasks. All containers should be pre-rinsed using G-RINSE™ in order to ensure there is no particulate matter or toxic substances residing.

### Supplementation of G-MOPS™, G-1™, G-2™, G-PGD™:

| Medium | HSA-solution™ | Final       |
|--------|---------------|-------------|
| [mL]   | volume [mL]   | volume [mL] |
| 9.5    | 0.5           | 10.0        |
| 19.0   | 1.0           | 20.0        |
| 28.5   | 1.5           | 30.0        |
| 38.0   | 2.0           | 40.0        |
| 47.5   | 2.5           | 50.0        |
| 57.0   | 3.0           | 60.0        |
| 66.5   | 3.5           | 70.0        |
| 76.0   | 4.0           | 80.0        |
| 85.5   | 4.5           | 90.0        |
| 95.0   | 5.0           | 100.0       |

### Supplementation of G-IVF™, G-FreezeKIt Blast™, G-ThawKIt Blast™:

| Medium | HSA-solution™ | Final       |
|--------|---------------|-------------|
| [mL]   | [mL]          | volume [mL] |
| 9.0    | 1.0           | 10.0        |
| 18.0   | 2.0           | 20.0        |
| 27.0   | 3.0           | 30.0        |
| 36.0   | 4.0           | 40.0        |
| 45.0   | 5.0           | 50.0        |
| 54.0   | 6.0           | 60.0        |
| 63.0   | 7.0           | 70.0        |
| 72.0   | 8.0           | 80.0        |
| 81.0   | 9.0           | 90.0        |
| 90.0   | 10.0          | 100.0       |

The cryopreservation solutions can be supplemented directly into the dishes. Simply place 900 µl of each cryopreservation solution into separate wells + 100 µl of G-MM™ or 100 µl HSA-solution™. It is recommended only to supplement the volume of medium expected to be used in one day.

### Specifications

|                                          |                     |
|------------------------------------------|---------------------|
| Sterile filtered                         | SAL 10 <sup>3</sup> |
| Mouse Embryo Assay (1-cell)              |                     |
| % expanded blastocysts (within 96 hours) | ≥ 80                |
| Bacterial endotoxins (LAL assay) [EU/mL] | < 0.5               |

LOT specific test results are available on the Certificate of Analysis provided with each delivery.

### Precautions

Discard product if bottle integrity is compromised. Do not use HSA-solution™ if it appears cloudy.

HSA-solution™ contains human serum albumin.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested for antibodies to HIV, HBc, HCV, and HTLV III and when not-reactive for HbsAg, HCV RNA and HIV-1 RNA and syphilis. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

Re-use may result in microbiological contamination and/or property changes in the product.

To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique.

The risk of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife's IVF media, have not been determined and are uncertain.

Not for injection.

Caution: Federal (US) law restricts this device to sale by or on behalf of a physician.

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