HYASE™-10X

US

EN: Indication for use

Medium for removal of cumulus cells.

Product Description
HYASE™ is a physiological salt buffer containing hyaluronidase, human serum albumin and gentamicin as an antibacterial agent.

For use after dilution 1:10 with either of the following two alternatives:

- G-MOPS™ PLUS or G-MOPS™ supplemented with G-MM™ or HSA-solution™ and equilibration at +37°C and ambient atmosphere.
- G-GAMETE™ and equilibration at +37°C and 6 % CO₂ atmosphere.

The concentration of hyaluronidase in HYASE is 800 IU/ml, before dilution.

Storage instructions and stability

Store dark at +2 to +8 °C.

HYASE is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

Directions for use

Dilute HYASE either with supplemented G-MOPS-G-MOPS PLUS/G-GAMETE to a ratio of 1:10, mix well.

Prepare and label sterile, non toxic multi-well dishes and prepare 3 wells of either supplemented G-MOPS-G-MOPS PLUS/G-GAMETE for every well of HYASE.

Pre-warm the multi well dishes with diluted HYASE and supplemented G-MOPS/G-MOPS PLUS for 15 minutes.

If G-GAMETE is used as diluting medium, equilibrate the dishes at +37°C and 6 % CO₂ for 3 hours.

Using a sterile, large bore pipette, place 3–5 oocytes into the HYASE. Gently pipette the HYASE and the oocytes to disperse the cumulus cells. Do not expose oocytes to HYASE solution for more than 30 seconds.

Use a sterile pipette to transfer oocytes from HYASE to the first washing well. Using a fine bore denudation pipette, aspirate each oocyte singly up and down to remove the corona. The diameter of the pipette should be slightly larger than that of the oocyte. Using a pipette that is too narrow or has a jagged edge may damage the oocyte.

Use the remaining two washing wells to rinse the oocyte. Transfer oocyte to new supplemented G-MOPS-G-MOPS PLUS/G-GAMETE or other appropriate medium for oocyte handling or, if needed, to culture medium for further incubation and maturation.

Specifications

Asepticly filtered	
Bacterial endotoxins (LAL assay) [EU/ml]	< 0.5*
*after diluting 1:10 in G-MOPS PLUS	
Cumulus Cell Removal Assay	PASS
Mouse Embryo Assay	Not performed
LOT specific test results are available on the Certificate of Analysis provided with each delivery.	

Precautions

Discard product if bottle integrity is compromised. Do not use HYASE if it appears cloudy.

HYASE contains hyaluronidase, human serum albumin and gentamicin. Do not use in patients with known hypersensitivity/allergy to any of the components.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested for antibodies to HIV, HBs, HCV, and HTLV III and non-reactive for HbsAg, HCV RNA and HIV-1 RNA and syphilis. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

Re-use may result in microbiological contamination and/or property changes in the product.

To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique. The risks of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife's IVF media, have not been determined and are uncertain.

Caution: Federal (US) law restricts this device to sale by or the order of a physician.

Description of ISO Symbols

STERILE LA	Sterilized using aseptic processing techniques	
→/←	Temperature limit	
Ⓢ	Do not re-use, discard after procedure.	
☞	Date of manufacture.	
📅	Use-by date.	
⚠	Caution: Consult accompanying documents.	
REF	Catalog number	
LOT	Batch code	
CE 2460	CE Mark (Conformité Européen)	

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Use the remaining two washing wells to

rinse the oocyte. Transfer oocyte to new supplemented G-MOPS-G-MOPS PLUS/G-GAMETE or other equivalent droplets for oocyte handling or, if needed, to culture medium for further incubation and maturation.

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To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique.

The risks of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife's IVF media, have not been determined and are uncertain.

Caution: Federal (US) law restricts this device to sale by or on the order of a physician.

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