## EN: Indication for use

Medium for gradient sperm separation.

## Product Description

All SpermGrad™ media are bicarbonate and HEPES buffered media containing silane-coated, colloid silica particles.

SpermGrad is available as:

-SpermGrad™ REF 10102 (30 mL) and 10099 (125 mL). Media for dilution into desired concentrations.

For use after dilution with either supplemented G-IVF™, G-IVF™ PLUS or with SpermRinse™.

For details on how to perform the protein supplementation see Package Inserts for G-MM™ or HSA-Solution™.

-SpermGrad™ REF 10138 (2x30 mL), 10238 (30 mL), 10338 (30 mL), 10139 (2x125 mL), 10239 (125 mL) and 10339 (125 mL) containing: Gradient solutions SpermGrad™ Upper Layer and/or SpermGrad™ Lower Layer. Ready to use solutions for gradient sperm separation.

SpermGrad™ is antibiotic-free while the media for dilution and sperm wash, such as supplemented G-IVF™, G-IVF™ PLUS or SpermRinse™, contain gentamicin as an antibacterial agent.

## Storage instructions and stability

Store dark at +2 to +8°C. SpermGrad™ is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

## Directions for use

Dilution: SpermGrad™ REF 10099 and 10102: Mix SpermGrad™ with supplemented G-IVF™, G-IVF™ PLUS or SpermRinse™ in separate tubes to obtain 90% and 45% stock solutions. For 90% stock solution, mix 9.0 mL SpermGrad™ with 1.0 mL supplemental G-IVF™, G-IVF™ PLUS or SpermRinse™.

SpermGrad™ REF 10138 (4.5 mL) with 5.5 mL SpermGrad™ Upper Layer and/or SpermGrad™ Lower Layer. Mix the solutions thoroughly. Store in sterile non-toxic tubes or sterile tissue culture flasks. Label with date and refrigerate until use.

SpermGrad™ REF 10138, 10238, 10338, 10139, 10239 and 10339: Ready to use.

Preparation of sample: Allow the solutions to reach ambient temperature before use.

Layer the gradient in a sterile, conical centrifuge tube marked with patient ID. Pipette 1.5 mL of 90 % solution (SpermGrad™ Lower Layer) into the tube first. Then add the slowest 1.5 mL of 45 % solution (SpermGrad™ Upper Layer) on top. Finally, the semen is gently layered on top of the gradient. Centrifuge the tube for 20 minutes at 300–600 g. Remove the two top layers. Transfer the sperm pellet with as little of the 90 % solution (SpermGrad™ Lower Layer) as possible to a sterile conical tube with 5 mL of equilibrated supplemented G-IVF™, G-IVF™ PLUS or SpermRinse™. Mix the solutions thoroughly. Store in sterile non-toxic tubes or sterile tissue culture flasks. Label with date and refrigerate until use.

SpermGrad™ REF 10138, 10238, 10338, 10139, 10239 and 10339: Ready to use.

Preparation of sample: Allow the solutions to reach ambient temperature before use.

Layer the gradient in a sterile, conical centrifuge tube marked with patient ID. Pipette 1.5 mL of 90 % solution (SpermGrad™ Lower Layer) into the tube first. Then add the slowest 1.5 mL of 45 % solution (SpermGrad™ Upper Layer) on top. Finally, the semen is gently layered on top of the gradient. Centrifuge the tube for 20 minutes at 300–600 g. Remove the two top layers. Transfer the sperm pellet with as little of the 90 % solution (SpermGrad™ Lower Layer) as possible to a sterile conical tube with 5 mL of equilibrated supplemented G-IVF™, G-IVF™ PLUS or SpermRinse™. Mix the solutions thoroughly. Store in sterile non-toxic tubes or sterile tissue culture flasks. Label with date and refrigerate until use.

SpermGrad™ REF 10138 (2x30 mL) and 10099 (125 mL). Media for dilution into desired concentrations.

For details on how to perform the protein supplementation see Package Inserts for G-MM™ or HSA-Solution™.

-SpermGrad™ REF 10102 (30 mL) and 10099 (125 mL). Human sperm recovery assay (recovery rate after gradient separation [%]) ≥ 25

Mouse Embryo Assay not performed

LOT-specific test results are available on the Certificate of Analysis provided with each delivery.

## Precautions

Discard product if bottle integrity is compromised. Do not use SpermGrad™ if it appears cloudy.

SpermGrad™ does not contain antibiotics.

Re-use may result in microbiological contamination and/or property changes in the product.

To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique.

The risks of reproductive toxicity and developmental toxicity for IVF media, including vitrolife's IVF media, have not been determined and are uncertain.

Not for injection.

Caution: Federal (US) law restricts this device to sale by or on the order of a physician.

## Description of ISO Symbols

STERILE[A] Sterilized using aseptic processing techniques

...K<sup>ac</sup> Temperature limit

(⊗) Do not re-use, discard after procedure.

...M<sup>W</sup> Date of manufacture

...U Use-by date

...A Caution: Consult accompanying documents.

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