EmbryoGlue®

US

en: Indication for use

Medium for embryo transfer.

Contraindications

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Medium for embryo transfer.

Contraindications

EmbryoGlue® contains recombinant human albumin (produced in yeast) and gentamicin. Do not use in patients with known hypersensitivity/allergy to any of the components.

Product Description

EmbryoGlue is a bicarbonate buffered medium containing recombinant human albumin, hyaluronan and gentamicin as an antibacterial agent.

For use after equilibration at +37°C and 6 % CO2 atmosphere.

The summary of safety and clinical performance can be found at [www.vitrolife.com](#)

Storage instructions and stability

Store dark at +2 to +8 °C.

EmbryoGlue is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

Directions for use

The product shall be used by an IVF professional.

The patient target group is an adult or reproductive-age population that undergoes fertility treatment.

Preparation of dishes

Add no less than 0.5mL of EmbryoGlue to the well and no less than 0.5 mL to the moat of a 1-well dish. Equilibrate the dish at +37°C and 6 % CO2 for up to 24 hours with the lid on.

In a non-humidified environment, dishes should not be stored for > 6 h.

Preparation of embryos

Equilibrate the embryos in EmbryoGlue for a minimum of 10 minutes at +37°C in a 6% CO2 environment prior to transfer. The embryos can stay in EmbryoGlue for up to 4 hours. To avoid evaporation, keep the lid on the dish.

In a non-humidified environment, complete the embryo transfer procedure within 6 h after dish preparation.

Preparation of transfer catheter

Attach the transfer catheter to a 1 mL, non-toxic syringe and draw up approximately 0.5 mL of the equilibrated medium from the moat. Rinse the catheter by expelling EmbryoGlue back into the moat.

Embryo transfer

After rinsing, draw approximately 0.1 mL of EmbryoGlue from the center well and expel into the moat until approximately 20 µL is left in the syringe. Under microscopic control, gently load the embryos into the catheter in approximately 5–10 µL of additional EmbryoGlue followed by a small amount of air. (The small pocket of air at the tip allows better imaging of the tip for ultrasound guided transfer.) For the embryo transfer, pass the tip of the catheter into the uterus approximately 1 cm from the top of the cavity and expel the embryos in a total volume of approximately 25–30 µL of medium. **Slowly withdraw the catheter while maintaining steady pressure on the plunger of the syringe.** Make a final microscopic examination of the catheter.

Specifications

Asepticlly filtered

Mouse Embryo Assay, 1-cell [% embryos developed to expanded blastocyst at 96 hours] ≥ 80

Bacterial endotoxins (LAL assay) [EU/mL] < 0.25

pH tested

Osmolality tested

LOT specific test results are available on the Certificate of Analysis provided with each delivery.

Precautions

Discard product if bottle integrity is compromised. Do not use EmbryoGlue if it appears cloudy.

Re-use may result in microbiological contamination and/or property changes in the product.

To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique.

The risk of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife’s IVF media, have not been determined and are uncertain.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Not for injection.

Discard product according to standard clinical practice for medical hazardous waste when the procedure is finished.

Caution: Federal (US) law restricts this device to sale by or on the order of a physician or practitioner trained in its use (Rx only).

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