## EmbryoSlide culture dishes

First issue 2022.11.03, revised 2023.07.07

### English

#### Scope of application

The EmbryoSlide culture dish is intended for the culture of oocytes and embryos in the EmbryoScope incubator.

#### Contraindications

None known.

#### Intended users

Embryologists, other laboratory personnel and clinic staff at IVF clinics trained by Vitrolife A/S-certified instructors.

#### Intended patient group

Patients undergoing infertility treatment.

#### Clinical benefits

- Improved embryo development
- Improved implantation/pregnancy rate
- Reduced pregnancy loss rate.

#### Limitations and restrictions on use

The EmbryoSlide culture dish is only intended for use in the EmbryoScope incubator and must be handled by trained personnel according to the instructions in this insert. Any serious incident that has occurred in relation to the culture dish should be reported to Vitrolife and to the competent authority of the Member State in which the user is established.

The EmbryoSlide culture dishes and lids are sterile. The pouches must only be opened in a sterile laminar flow hood, and the culture dish must be covered by the lid when in use.

#### Warnings and precautions

##### WARNINGS

The EmbryoSlide culture dish is intended for single use only and MAY NOT be re-used. Any attempt by the user to clean and re-sterilise the culture dish may result in contamination with microorganisms or other risks of device failure.

To avoid contamination with microorganisms, always place the culture dish in a sterile laminar flow hood while loading and generally handling the culture dish.

Always label the culture dish appropriately and validate the label when the embryo is transferred either to a new device or to the patient.

#### PRECAUTIONS

Do not use the rinsing wells for incubating embryos as no images are acquired from these wells.

The reservoir and rinsing wells MUST always be covered by a confluent oil layer of 1.4 ml of IVF-grade oil. This prevents the medium from evaporating during incubation.

Evaporation of the medium can change osmolality, which may affect embryo development.

Always place the culture dish in a sterile laminar flow hood while loading it. Take care while loading and handling the culture dish to reduce the risk of spilling any oil or medium.

If ANY oil or medium is accidentally spilled from the culture dish while it is loaded or handled, the embryos MUST be transferred to a new culture dish to prevent them from being adversely affected by the spillage.

In case any oil is spilled on the culture dish outside of the reservoir, gently wipe off the oil by using a lint-free tissue paper. If any oil is present on the lid, the lid must be replaced.

There is a slight variation in how much the temperature decreases in the microwells and rinsing wells during handling:

- First two minutes: The temperature decreases less than 0.1°C in the microwells and approx. 0.7°C in the rinsing wells.
- Between 2 and 5 minutes: The temperature decreases approx. 0.2°C in the microwells and 0.6°C in the rinsing wells. Always verify the temperature on your own heating stage.

If any bubbles are present in the well or oil layer after preparation, they may block the camera light and compromise image quality. In rare cases, bubbles may dislodge the embryo from the bottom of the microwell.

Any bubbles must be removed carefully and relatively quickly to avoid evaporation of medium.

Take care not to remove any medium when removing the bubbles.

Discard the product according to standard clinical practice for medical hazardous waste when the procedure is finished.

#### Product description

The EmbryoSlide culture dish is a single-use sterile dish designed for the culture of embryos in the EmbryoScope incubator.

The EmbryoSlide culture dish contains two types of wells: 12 culture wells in which the embryos reside during incubation and four rinsing wells that are used for rinsing and handling the embryos. The embryos reside in a central depression –

the microwell – inside each culture well. The culture wells are cylindrical and have a diameter of 4.09 mm. The microwells are conical with a flat bottom and have a bottom diameter of 0.25 mm. See the illustration on the last page of this insert for an overview of the culture dish.

#### Preparing the culture dish

Prepare one culture dish at a time to minimise handling time. Work on a non-heated workbench and use cold medium and oil.

1. Remove the culture dish from the pouch in a sterile laminar flow hood.
2. Fill all microwells with medium using a micropipette (max. diameter: 200 µm). One filling of the tip will suffice to completely fill all microwells. The tip of the pipette must touch the side of the microwell during filling. Slightly overfill each microwell.
3. Use a standard pipette to immediately fill an additional 25 µl of medium into each well. It is important to release the medium from the pipette close to the bottom of the microwell at a constant, not too fast rate.
4. Fill each rinsing well with 25 µl of medium.
5. Load 1.4 ml of oil certified for IVF into the reservoir. Make sure that all culture wells and rinsing wells are covered with a confluent oil layer. Add an additional amount of oil for each well or reservoir not filled with medium.
6. Push up any large bubbles with a micropipette and then remove them.
7. Cover the culture dish with the lid and equilibrate it overnight.
8. Identify and remove any bubbles under a stereo microscope.
9. Load embryos into the centre of the microwells using a micropipette.
10. Place the culture dish in the EmbryoScope incubator.

#### Product name

Culture dishes.

#### Product model and specification

Model and specification: FT-S-ES-D, EmbryoSlide.

#### Product registration certificate number/product technical requirement number

Product registration certificate number:

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Product technical requirement number:

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#### Product features

|   |                      |
|---|----------------------|
| Sterile   | SAL 10 <sup>-6</sup> |
| Bacterial endotoxin (LAL experiment) [EU/item]  | < 20                 |
| Mouse in-vitro embryo test (1-cytoplasm mouse embryo) [expanded blastocyst percentage on the fifth day] | ≥ 80%                |

All products include the quality inspection analysis report for the corresponding batch number.

#### Product main structure and composition

This product is made of uncoated polystyrene.

#### Manufacturing date, expiry date and storage and transportation conditions

Manufacturing date: see the product label.  
Expiry date: see the product label.  
Storage conditions: store at room temperature.  
Transportation conditions: transport at room temperature.

#### Registration applicant/manufacturing organisation, domicile, manufacturing address, production licence number and contact information

Name: Vitrolife A/S  
Domicile and manufacturing address: Jens Juuls Vej 20, 8260 Viby J, Denmark  
Production licence number: CVR (DK) 27406793  
Telephone: +45 7221 7900  
Fax: +45 7221 7901

#### Agent/after-sales service company name, address and contact information

Name: Vitrolife (Beijing) Medical Device Co., Ltd.  
Address: Room 2006, 20th (17) Floor, Building No. 3, No. 6 Futong East Street, Chaoyang District, Beijing, China  
Zip code: 100102  
Telephone: +86 10 64036613  
Fax: +86 10 64036613  
E-mail: support.asia@vitrolife.com  
www.vitrolife.com

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

E-mail: support.asia@vitrolife.com

www.vitrolife.com

| Symbol  | Title  | Description   | 描述                                     |
|---|--|---|--|
|    | <b>Catalogue number</b><br>ISO 15223-1:2021<br>5.1.6                   | Indicates the manufacturer's catalogue number so that the medical device can be identified.                                       | 表示制造商的目录编号、以便识别医疗设备。                   |
|    | <b>Batch code</b><br>ISO 15223-1:2021<br>5.1.5                         | Indicates the manufacturer's batch code so that the batch or lot can be identified.   | 表示制造商的批次代码、以便识别批次。                     |
|    | <b>Use by date</b><br>ISO 15223-1:2021<br>5.1.4                        | Indicates the date after which the medical device is not to be used.  | 表示不能再使用该医疗设备的日期。                       |
|    | <b>Do not re-use</b><br>ISO 15223-1:2021<br>5.4.2                      | Indicates a medical device that is intended for one use or for use on a single patient during a single procedure.                 | 表示某个医疗设备为一次性使用、或在一个手术中单独用于某个患者。        |
|    | <b>Do not use if packaging is damaged</b><br>ISO 15223-1:2021<br>5.2.8 | Indicates a device that should not be used if the package has been damaged or opened.   | 表示某个医疗设备如果包装损坏或打开、则不能使用。               |
|    | <b>CE mark</b>   | Declaration by the manufacturer that the device meets all applicable requirements in the Medical Device Regulation (EU) 2017/745. | 制造商声明该设备符合医疗器械法规（欧盟）2017/745 中的所有适用要求。 |
|    | <b>Sterilised using irradiation</b><br>ISO 15223-1:2021<br>5.2.4       | Indicates a medical device that has been sterilised using irradiation.  | 表示医疗设备已辐射消毒。                           |
|    | <b>Consult the instructions for use</b><br>ISO 15223-1:2021<br>5.4.3   | Indicates the need for the user to consult the instructions for use.  | 表示用户需要查阅使用说明。                          |
|    | <b>Temperature limit</b><br>ISO 15223-1:2021<br>5.3.7                  | Indicates the temperature limits to which the medical device can be safely exposed.   | 表示医疗设备可以安全暴露的温度范围。                     |
|    | <b>Manufacturer</b><br>ISO 15223-1:2021<br>5.1.1                       | Indicates the medical device manufacturer.  | 表示医疗设备的制造商。                            |
|   | <b>Date of manufacture</b><br>ISO 15223-1:2021<br>5.1.3                | Indicates the date when the medical device was manufactured.  | 表示医疗设备的制造日期。                           |
|  | <b>Medical device</b><br>ISO 15223-1:2021<br>5.7.7                     | Indicates that the device is a medical device.  | 表示该设备是医疗器械。                            |
|  | <b>Unique device identifier</b><br>ISO 15223-1:2021<br>5.7.10          | Indicates a carrier that contains unique device identifier information.   | 表示包含唯一设备标识信息的载体。                       |

| No.     | 1         | 2         | 3         | 4            | 5    | 6            | 7      | 8         | 9    |
|---------|-----------|-----------|-----------|--------------|------|--------------|--------|-----------|------|
| English | Reservoir | Microwell | Label tab | Rinsing well | Well | Rinsing well | Handle | Microwell | Well |
| 中文      | 储油孔       | 微孔        | 标签        | 冲洗孔          | 培养孔  | 冲洗孔          | 手柄     | 微孔        | 培养孔  |

**Technical support**

Tel.: +45 7023 0500  
 E-mail: support.embryoscope@vitrolife.com  
 www.vitrolife.com

 Vitrolife A/S  
 Jens Juuls Vej 20  
 DK-8260 Viby J  
 Tel.: +45 7221 7900