EmbryoSlide Flex culture dishes

First issue 2019.04.17, revised 2022.03.24

EN (UK)

Intended use

The EmbryoSlide Flex culture dish is an accessory for embryo storage related to the EmbryoScope Flex incubator, which provides an environment with controlled temperature and gas concentrations (CO₂ and optionally O₂) for the development of gametes and/or embryos at or near body temperature.

Indications for use

The EmbryoSlide Flex culture dish is intended for preparing, storing and transferring human gametes and/or embryos. The EmbryoSlide Flex culture dish must be used together with the EmbryoScope Flex incubator.

Contraindications

None known.

Intended users

Embryologists, other laboratory personnel and clinic staff at IVF clinics trained by Vitrolife A/S-certified instructors.

Intended patient group

Patients undergoing infertility treatment.

Clinical benefits

- Improved embryo development
- Improved implantation/pregnancy rate
- Reduced pregnancy loss rate.

Limitations and restrictions on use

The EmbryoSlide Flex culture dish is only intended for use in the EmbryoScope Flex incubator and must be handled by trained personnel according to the instructions in this insert. Any serious incident that has occurred in relation to the culture dish should be reported to Vitrolife and to the competent authority of the Member State in which the user is established.

The EmbryoSlide Flex culture dishes and lids are sterile. The pouches must only be opened in a sterile laminar flow hood, and the culture dish must be covered by the lid when in use.

Warnings and precautions

WARNINGS

The EmbryoSlide Flex culture dish is intended for single use only and MAY NOT be re-used. Any attempt by the user to clean and re-sterilise the culture dish may result in contamination with microorganisms or other risks of device failure.

To avoid contamination with microorganisms, always place the culture dish in a sterile laminar flow hood while loading and generally handling the culture dish.

Always label the culture dish appropriately and validate the label when the embryo is transferred either to a new device or to the patient.

PRECAUTIONS

Do not use the rinsing wells for incubating embryos as no images are acquired from these wells.

The reservoir and rinsing wells MUST always be covered by a confluent oil layer of 1.2 ml of IVF-grade oil. This prevents the medium from evaporating during incubation.

Evaporation of the medium can change osmolality, which may affect embryo development.

Always place the culture dish in a sterile laminar flow hood while loading it. Take care while loading and handling the culture dish to reduce the risk of spilling any oil or medium.

If ANY oil or medium is accidentally spilled from the culture dish while it is loaded or handled, the embryos MUST be transferred to a new culture dish to prevent them from being adversely affected by the spillage.

In case any oil is spilled on the culture dish outside of the reservoirs, gently wipe off the oil by using a lint-free tissue paper. If any oil is present on the lid, the lid must be replaced.

If any bubbles are present in the well or oil layer after preparation, they may block the camera light and compromise image quality. In rare cases, bubbles may dislodge the embryo from the bottom of the microwell.

Any bubbles must be removed carefully and relatively quickly to avoid evaporation of medium.

Take care not to remove any medium when removing the bubbles.

Discard the product according to standard clinical practice for medical hazardous waste when the procedure is finished.

Product description

The EmbryoSlide Flex culture dish is a single-use sterile dish designed for the culture of embryos in the EmbryoScope Flex incubator.

The EmbryoSlide Flex culture dish contains two types of wells: six culture wells in which the embryos reside during incubation and two rinsing wells that are used for rinsing and handling the embryos. The embryos reside in a central depression – the microwell – inside each culture well. The culture wells are cylindrical with a diameter of 1.5 mm and are placed in a common medium reservoir. The microwells are conic with a curved bottom and have a bottom diameter of 0.28 mm. See the illustration on the last page of this insert for an overview of the culture dish.

Preparing the culture dish

Prepare one culture dish at a time to minimise handling time. Work on a non-heated workbench and use cold medium and oil.

1. Remove the culture dish from the pouch in a sterile laminar flow hood.
2. Fill all microwells with medium using a micropipette (max. diameter: 200 µm). One filling of the tip will suffice to completely fill all microwells. To avoid bubbles: Slightly overfill each microwell to avoid creating a depression in the surface of the medium. The tip of the pipette must touch the side of the microwell during filling.
3. Use a standard pipette to immediately fill an additional 140 µl of medium into the wells in the reservoir. Carefully slide the tip over the wells while releasing the medium. Make sure to fill the reservoir completely.
4. Fill each rinsing well with 35 µl of medium.
5. Load 1.2 ml of oil certified for IVF into the reservoir. Make sure that all culture wells and rinsing wells are covered with a confluent oil layer. Add an additional amount of oil for each well or reservoir not filled with medium.
6. Push up any large bubbles with a micropipette and then remove them.
7. Cover the culture dish with the lid and equilibrate it overnight.
8. Identify and remove any bubbles under a stereo microscope.
9. Place the barcode label on the dedicated label tab on the culture dish. The barcode must be placed so that it is possible to read the text on the label when the culture dish is correctly inserted into the incubator. Make sure that the label is smoothed out as best as possible to avoid wrinkles in the paper. If the barcode is in any way damaged or not smoothly applied, the label cannot be read. In this case, print a new label from the EmbryoViewer software.
10. Load embryos into the centre of the microwells using a micropipette.
11. Place the culture dish in the EmbryoScope Flex incubator.

EN (US)

Intended use

The EmbryoSlide Flex culture dish is an accessory for embryo storage related to the EmbryoScope Flex incubator, which provides an environment with controlled temperature and CO₂ for the development of embryos at or near body temperature.

Indications for use

The EmbryoSlide Flex culture dish is intended for preparing, storing and transferring human embryos. The EmbryoSlide Flex culture dish must be used together with the EmbryoScope Flex incubator.

Restrictions on use

The EmbryoSlide Flex culture dish is only intended for use in the EmbryoScope Flex incubator and must be handled by trained personnel according to the instructions in this insert.

The EmbryoSlide Flex culture dishes and lids are sterile. The pouches must only be opened in a sterile laminar flow hood, and the culture dish must be covered by the lid when in use.

Warnings and precautions

WARNINGS

The EmbryoSlide Flex culture dish is intended for single use only and MAY NOT be re-used. Any attempt by the user to clean and re-sterilize the culture dish may result in contamination with microorganisms or other risks of device failure.

To avoid contamination with microorganisms, always place the culture dish in a sterile laminar flow hood while loading and generally handling the culture dish.

Always label the culture dish appropriately and validate the label when the embryo is transferred either to a new device or to the patient.

PRECAUTIONS

Federal law restricts this device to sale by or on the order of a physician or a practitioner trained in its use.

Do not use the rinsing wells for incubating embryos as no images are acquired from these wells.

The reservoir and rinsing wells MUST always be covered by a confluent oil layer of 1.2 ml of IVF-grade oil. This prevents the medium from evaporating during incubation.

Evaporation of the medium can change osmolality, which may affect embryo development.

Always place the culture dish in a sterile laminar flow hood while loading it to keep both hands free and reduce the risk of spilling any oil or medium.

If ANY oil or medium is accidentally spilled from the culture dish while it is loaded or handled, the embryos MUST be transferred to a new culture dish to prevent them from being adversely affected by the spillage.

In case any oil is spilled on the culture dish outside of the reservoirs, gently wipe off the oil by using a lint-free tissue paper. If any oil is present on the lid, the lid must be replaced.

If any bubbles are present in the well or oil layer after preparation, they may block the camera light and compromise image quality. In rare cases, bubbles may dislodge the embryo from the bottom of the microwell.

Any bubbles must be removed carefully and relatively quickly to avoid evaporation of medium.

Take care not to remove any medium when removing the bubbles.

Preparing the culture dish

Prepare one culture dish at a time to minimize handling time. Work on a non-heated workbench and use cold medium and oil.

1. Remove the culture dish from the pouch in a sterile laminar flow hood.
2. Fill all microwells with medium using a micropipette (max. diameter: 200 µm). One filling of the tip will suffice to completely fill all microwells. To avoid bubbles: Slightly overfill each microwell to avoid creating a depression in the surface of the medium. The tip of the pipette must touch the side of the microwell during filling.
3. Use a standard pipette to immediately fill an additional 140 µl of medium into the wells in the reservoir. Carefully slide the tip over the wells while releasing the medium. Make sure to fill the reservoir completely.
4. Fill each rinsing well with 35 µl of medium.
5. Load 1.2 ml of oil certified for IVF into the reservoir. Make sure that all culture wells and rinsing wells are covered with a confluent oil layer. Add an additional amount of oil for each well or reservoir not filled with medium.

6. Push up any large bubbles with a micropipette and then remove them.
7. Cover the culture dish with the lid and equilibrate it overnight.
8. Identify and remove any bubbles under a stereo microscope.
9. Place the barcode label on the dedicated label tab on the culture dish. The barcode must be placed so that it is possible to read the text on the label when the culture dish is correctly inserted into the incubator. Make sure that the label is smoothed out as best as possible to avoid wrinkles in the paper. If the barcode is in any way damaged or not smoothly applied, the label cannot be read. In this case, print a new label from the EmbryoViewer software.
10. Load embryos into the centre of the microwells using a micropipette.
11. Place the culture dish in the EmbryoScope Flex incubator.

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