G-PGD™

US

en: Indication for use

Medium for embryo biopsy.

Contraindications

G-PGD™ contains gentamicin. Do not use in patients with known hypersensitivity/allergy to the component.

Product Description

G-PGD is a calcium and magnesium free MOPS buffered medium containing gentamicin as an antibacterial agent. For use after the addition of HSA-solution™ and equilibration at +37 °C and ambient atmosphere.

The summary of safety and clinical performance can be found at www.vitrolife.com

Storage instructions and stability

Store dark at +2 to +8 °C.

G-PGD is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

For Non-UE: Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

EU: Media bottles can be used for up to two weeks after first opening, use aseptic technique and minimize the time outside the refrigerator. Record the opening date on the bottle. Discard excess media no later than two weeks after first opening.

Directions for use

The product shall be used by an IVF professional. The patient target group is an adult population that undergoes IVF treatment including PGT of embryos.

Media to be supplemented should be aliquoted into sterile non-toxic tissue culture grade tubes or flasks. For details on how to perform the protein supplementation see Package Inserts for HSA-solution.

G-PGD must never be placed in the CO₂ incubator but rather warmed in an incubator without gas or in a heating block. If oil overlay is used, the oil should not be equilibrated in CO₂.

Embryo biopsy is usually performed on the morning of day 3, before compaction of the blastomeres begins. A hole in the zona can be drilled either with Tyrodé's acid or with a laser. A separate pipette is needed to remove the biopsied blastomere.

Prepare dishes with droplets of pre-warmed supplemented G-PGD covered with OVOIL™.

Place the embryo in supplemented G-PGD at +37 °C. Secure the embryo with the holding pipette.

Drill a hole in the zona. Introduce a clean biopsy pipette into the hole and gently aspirate one or two blastomeres for analysis.

The biopsied embryo should be placed into equilibrated G-2™ PLUS for culture until the results of the embryo biopsy analysis are confirmed. The non-affected embryos may be selected for fresh embryo transfer on day 5 or cryopreservation if there is a surplus. The embryo transfer is usually performed on day 5, after the genetic diagnosis has been completed.

Caution: Spend time adjusting the embryo position on the holding pipette so that a nucleated blastomere can be clearly identified and biopsied. If compaction has occurred, there will be junctions between the blastomeres and they will be more difficult to remove. Be patient, pull and move the pipette across and back to ease the blastomere out.

Specifications

Asepticly filtered

Mouse Embryo Assay, 1-cell [% embryos developed to expanded blastocyst at 96 hours] ≥ 80

Bacterial endotoxins (LAL assay) [EU/mL] < 0.25

pH tested

Osmolality tested

LOT specific test results are available on the Certificate of Analysis provided with each delivery.

Precautions

Discard product if bottle integrity is compromised. Do not use G-PGD if it appears cloudy.

To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique.

The risk of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife's IVF media, have not been determined and are uncertain.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Not for injection.

Discard the product according to standard clinical practice for medical hazardous waste when the procedure is finished.

Caution: Federal (US) law restricts this device to sale by or on the order of a physician or practitioner trained in its use (Rx only).

Symbols glossary

* Medical devices – Symbols to be used with information to be supplied by the manufacturer (ISO 15223-1:2021).

Symbol	Title	Reference and symbol ID
	Consult instructions for use	*Ref no: 5.4.3
	Caution	*Ref no: 5.4.4
	Manufacturer	*Ref no: 5.1.1
	Date of manufacture	*Ref no: 5.1.3
	Use-by date	*Ref no: 5.1.4
	Batch code	*Ref no: 5.1.5
	Catalogue number	*Ref no: 5.1.6
	Sterilized using aseptic processing techniques	*Ref no: 5.2.2
	Single sterile barrier system	*Ref no: 5.2.1
	Keep away from sunlight (light sources)	*Ref no: 5.3.2
	Temperature limit	*Ref no: 5.3.7
	Contains a medicinal substance	*Ref no: 5.4.7
	Medical Device	*Ref no: 5.7.7
	European conformity (CE) mark with Notified Body identification number as defined in Regulation (EU) 2017/745 Annex V	

Caution: Federal (US) law restricts this device to sale by or on the order of a physician or practitioner trained in its use (Rx only).

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