HASE™-10X

US

0N: Indication for use

Medium for removal of cumulus cells.

Contraindications

HAASE-10X contains gentamicin. Do not use in patients with known hypersensitivity/allergy to the component.

Product Description

HAASE-10X is a physiological salt buffer containing hyaluronidase, human serum albumin and gentamicin as an antibacterial agent.

For use after dilution 1:10 with either of the following two alternatives:

1) G-MOPStm PLUS or G-MOPStm PLUS+ supplemented with HSA-solution** and equilibration at +37°C and ambient atmosphere.

2) G-GAMETEm and equilibration at +37°C and 6 % CO2 atmosphere.

The concentration of hyaluronidase in HAASE-10X is 800 IU/mL, before dilution.

The summary of safety and clinical performance can be found at www.vitrolife.com

Storage instructions and stability

Store dark +2 to +8 °C.

HAASE-10X is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

Directions for use

The product shall be used by an IVF professional. The patient target group is an adult or reproductive-age population that undergoes IVF treatment or fertility preservation.

Dilute HAASE-10X either with supplemented G-MOPStm-MOPS PLUS+G-GAMETE to a ratio of 1:10, mix well.

Prepare and label sterile, non-toxic multi-well dishes and prepare 3 wells of either supplemented G-MOPStm-MOPS PLUS+G-GAMETE for every well of HAASE. Pre-warm the multi well dishes with diluted HAASE and supplemented G-MOPStm-MOPS PLUS for 15 minutes. If G-GAMETE is used as diluting medium, equilibrate the dishes at +37 °C and 6 % CO2 for 3 hours.

Use a sterile pipette to transfer oocytes from HAASE to the first washing well. Gently pipette the HAASE and the oocytes to disperse the cumulus cells. Do not expose oocytes to HAASE solution for more than 30 seconds.

Use a sterile pipette to transfer oocytes from HAASE to the first washing well. Gently pipette the HAASE and the oocytes, aspirate each oocyte singly up and down to remove the corona. The diameter of the pipette should be slightly larger than that of the oocyte. Using a pipette that is too narrow or has a jagged edge may damage the oocyte.

Use the remaining two washing wells to rinse the oocyte. Transfer oocyte to new supplemented G-MOPStm-MOPS PLUS+G-GAMETE or other equivalent droplets for oocyte handling or, if needed, to culture medium for further incubation and maturation.

Specifications

Asepticly filtered Bacterial endotoxins (LAL assay) [EU/mL] Cumulus Cell Removal Assay Musem Embryo Assay Not performed pH tested Osmolality tested **after diluting 1:10 in G-MOPS PLUS	<0.5* PASS Not performed
---	---

LOT specific test results are available on the Certificate of Analysis provided with each delivery.

Precautions

Discard product if bottle integrity is compromised. Do not use HAASE-10X if it appears cloudy.

HAASE-10X contains hyaluronidase and human serum albumin.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested for antibodies to HIV, HBV, HCV, and HTLV III and non-reactive for HbsAg, HCV RNA and HIV-1 RNA and syphilis. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

Re-use may result in microbiological contamination and/or property changes in the product.

To avoid contamination Vitrolife strongly recommends that media should be opened and used with aseptic technique.

The risks of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife's IVF media, have not been determined and are uncertain.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the user and/or patient is established.

Not for injection.

Discard the product according to standard clinical practice for medical hazardous waste when the procedure is finished.

Symbols Glossary

** Medical devices – Symbols to be used with information to be supplied by the manufacturer (ISO 15223-1:2021).

Symbol	Title	Reference and symbol ID
	Consult instructions for use	**Ref no. 5.4.3
	Caution	**Ref no. 5.4.4
	Manufacturer	**Ref no. 5.1.1
	Date of manufacture	**Ref no. 5.1.3
	Use by date	**Ref no. 5.1.4
	Batch code	**Ref no. 5.1.5
	Catalogue number	**Ref no. 5.1.6
	Sterilized using aseptic processing techniques	**Ref no. 5.2.2
	Single sterile barrier system	**Ref no. 5.2.11
	Do not re-use	**Ref no. 5.3.2
	Temperature limit	**Ref no. 5.4.7
	Keep away from sunlight (light sources)	**Ref no. 5.3.2
	Contains a medicinal substance	**Ref no. 5.4.7
	Contains human blood or plasma derivatives	**Ref no. 5.4.6
	Contains biological material of animal origin	**Ref no. 5.4.8
	Medical Device	**Ref no. 5.7.7
	CE Mark (Conformité Européenne)	European conformity (CE) mark with Notified Body identification as defined in Annex IV of the EU 2017/745 Directive

Caution: (Federal US) law restricts this device to sale by or on the order of a physician or practitioner trained in its use (Rx only).

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