# Micromanipulation Pipettes

Vitrolife

US

## en: Indication for use

Micromanipulation pipettes are intended for micro-manipulation of embryos, cells and gametes during assisted reproductive technology (ART).

Biopsy Pipette - For aspiration of cells for pre-implantation genetic testing.

Handling Pipette - For handling and manipulation of embryos, cells and gametes during ART procedures.

Hatching Pipette - To create a hole in the zona pellucida to enable assisted hatching or embryo biopsy.

Holding Pipette - To hold the oocyte or embryo when performing ICSI or other micromanipulation procedures.

ICSI Pipette - For aspiration and intracytoplasmic injection of a single sperm into the oocyte.

PZD Pipette - To create a small slit in the zona pellucida to enable assisted hatching or embryo biopsy.

## Contraindications

None known.

## Product Description

Micromanipulation pipettes are sterile, disposable, borosilicate glass lab products. Handling Pipettes are equipped with a bulb.

## Storage Instructions

Store at +2 to +30°C.

## Directions for use

The product shall be used by an IVF professional.

The patient target group is an adult or reproductive-age population that undergoes fertility treatment.

Check the packaging for signs of any damage.

Peel open the lid on the pipette holder. Remove the lid completely. Grasp the pipette in the middle and pull the pipette straight up.

## Specifications

Sterilized using irradiation SAL 10<sup>6</sup>

Mouse Embryo Assay, 1-cell [% embryos developed to expanded blastocyst at 96 hours] ≥ 80

Bacterial Endotoxins (LAL assay) < 1.0 EU/device

LOT specific test results are available on the Certificate of Analysis provided with each delivery.

## Precautions

Discard product if the product or packaging is damaged or wet or if the pipette has been contaminated during unpacking or mounting. Do not use broken pipettes.

Glass-shards may be present in spite of the 100% visual inspection, always inspect pipette before use.

Micromanipulation pipettes are intended for single use only and SHALL NOT BE REUSED or RESTERILIZED. Reuse may cause contamination, poor embryo survival, and failed procedure.

Discard product according to standard clinical practice for medical hazardous waste and glass sharps when the procedure is finished.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## Symbols glossary

\*\* Medical devices – Symbols to be used with information to be supplied by the manufacturer (ISO 15223-1:2021).

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