RapidWarm Blast

US

en: Indication for use

Media for warming of vitrified human blastocyst embryos.

Contraindications

RapidWarm® Blast contains media for the warming of vitrified human blastocyst stage embryos. The media consists of a MOPS buffered solution containing gentamicin as an antibacterial agent and human serum albumin.

Product Description

SUPPLEMENTED WITH HSA

RapidWarm Blast contains three media for the warming of vitrified human blastocyst stage embryos. The media consists of a MOPS buffered solution containing gentamicin as an antibacterial agent and human serum albumin.

Storage instructions and stability

Store dark at +2 to +8 °C.

RapidWarm Blast is stable until the expiry date shown on the bottle labels and the LOT specific certificate of analysis.

Media bottles should not be stored after opening.

Discard excess media after completion of procedure.

Use of storage device

Use a legally marketed storage device indicated for use in human blastocyst vitrification procedures, in order to assure sufficient cooling and warming rates.

Use a closed storage system to prevent the potential risk of virus or other contamination of samples.

Perform the actual vitrification and warming according to the Instructions for Use for the storage device.

Directions for use

The product shall be used by an IVF professional.

The patient target group is an adult or reproductive-age patient who undergoes IVF treatment or fertility preservation.

The following is the general procedure for using RapidWarm Blast:

Note that the timing has to be accurate.

Warning:

Phase 0.5-1 ml of each of the following media into separate wells of a 4-well plate and warm to 37 °C.

Warm 1 Blast

Warm 2 Blast

Warm 3 Blast

All manipulations of the blastocysts are carried out at 37 °C (on the heated stage).

Remove the cryovial containing the vitrified blastocysts from the storage container. Follow the warming instructions for the specific cryovial used.

Immediately after warming, transfer blastocysts into Warm 1 Blast.

Add the blastocysts to fall from the device and sink to the bottom. Leave for 2 min.

Transfer the blastocysts into Warm 2 Blast and let the blastocysts remain in the solution for 5-10 min.

Transfer the blastocysts into Warm 3 Blast and let the blastocysts remain in the solution for 5-10 min.

Remove the blastocysts in culture media several times and allow them to acclimatise to laboratory practice. Follow the recommendations using EmbryoGlu® for this procedure.

Specifications

Media Embryo Assay, 1-cell [%] (re-expanded blastocysts 24 hours post-thaw) ≥ 70

Bacterial endotoxins (LAL assay) [EU/mL] ≤ 0.5 pH tested

Osmolality tested

LOT-specific test results are available on the Certificate of Analysis provided with each delivery.

Precautions

The long-term safety of vitrification and/or blastocyst collection on children born following this method of embryo cryopreservation procedure has not been established.

The risk of reproductive toxicity and developmental toxicity for IVF media, including VitroLife's IVF media, have not been determined and are uncertain.

The safety and effectiveness of vitrification have not been fully evaluated in human embryos that have not yet reached the blastocyst stage of development.

Discard product if bottle integrity is compromised. Do not use RapidWarm Blast if the product is compromised.

Discard product after human serum albumin is compromised.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested for antibodies to HIV, HBV, HCV, and HTLV 1/II and non-reactive for HbsAg, HCV RNA and HIV-1 RNA and syphilis. No known test methods can assure that products derived from human blood will not transmit infectious agents.

Re-use may result in microbiological contamination and/or property changes in the product.

To avoid contamination VitroLife strongly recommends that product should be opened and used only with aseptic technique.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Not for injection.

Discard the product according to standard clinical practice for sharps and contaminated equipment when the procedure is finished.

Symbol's glossary

* Medical devices – Symbols to be used with information to be supplied by the manufacturer (ISO 15223-1:2021).

Symbol	Title	Reference and symbol ID
	Consult instructions for use	*Ref no: 5.4.3
	Caution	*Ref no: 5.4.4
	Manufacturer	*Ref no: 5.4.1
	Date of manufacture	*Ref no: 5.1.3
	Use by date	*Ref no: 5.1.4
	Batch code	*Ref no: 5.1.5
	Catalogue number	*Ref no: 5.1.6
	Sterilized using aseptic processing techniques	*Ref no: 5.2.2
	Single sterile barrier system	*Ref no: 5.2.11
	Temperature limit	*Ref no: 5.3.3
	Keep away from sunlight (light sources)	*Ref no: 5.3.7
	Contains a medicinal substance	*Ref no: 5.4.7
	Contains human blood or plasma derivatives	*Ref no: 5.4.6
	Medical Device	*Ref no: 5.4.7.1
	European conformity (CE) mark with Notified Body identification number as defined in Regulation (EU) 2017/745 Annex V	*Ref no: 5.4.7.2

Caution: Federal (US) law restricts this device to sale by or on the order of a physician or practitioner trained in its use (R only).

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