

RapidWarm™ Omni contains gentamicin. Do not use in patients with known hypersensitivity/allergy to the component.

Product Description

SUPPLEMENTED WITH HSA
RapidWarm™ Omni contains four solutions for the warming of unfertilized oocytes through to blastocyst stage embryos. The solutions consist of a MOPS-based medium containing gentamicin as an antibacterial agent, and human serum albumin.

Warm™ 1 Omni contains sucrose
Warm™ 2 Omni contains sucrose
Warm™ 3 Omni contains sucrose
Warm™ 4 Omni contains no sucrose

For use after warming to +37 °C in ambient atmosphere. The summary of safety and clinical performance can be found on the product label.

Storage Instructions and stability

Storage date is at 2 to +8 °C.
RapidWarm™ Omni is stable until the expiry date shown on the container and the L-OTC-specific Certificate of Analysis.

Media bottles should be used for up to two weeks after first opening, use aseptic technique and minimize the time outside the refrigerator. Record opening date on the bottle. Discard excess media no later than 10 weeks after opening.

US: Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

Directions for use

The product is used by an IVF professional. The patient or user is to be informed of an reproductive population that undergoes IVF treatment or fertility preservation.

The following is the general procedure for using RapidWarm™ Omni.

Note: Timing with warming is critical, ensure you follow the protocol precisely.

Warning: Please use one of each of the following media into separate wells of a multi-well plate and warm to 37 °C in ambient atmosphere.

• Warm 1 Omni
• Warm 2 Omni
• Warm 3 Omni
• Warm 4 Omni

Note: The recommended volumes should not be changed. Volume changes will affect temperature control in the first warming solution as well as osmolality, which may result in damage to the embryo.

It is very important to keep the temperature of all warming solutions at 37 °C at all times. If more than one cryoprotectant is to be warmed in the same dish, make sure that the temperature reaches 37 °C after each warming procedure.

All manipulations of the oocytes or embryos are done at 37 °C (on the heated stage) in ambient atmosphere.

Remove the cryoprotectant containing the unfertilized oocytes or embryos from the cryotank or the warming stage and place the specific cryoprotectant.

Immediately, transfer oocytes or embryos into Warm 1 Omni. It is important to use 1 ml of Warm 1 Omni to ensure the temperature of the medium is not perturbed by the addition of the cryoprotectant.

Allow the oocytes or embryos to fall from the device and sink to the bottom. Leave for 1 minute.

Transfer the oocytes or embryos into Warm 2 Omni and let the oocytes or embryos remain in the solution for 3 minutes.

Transfer the oocytes or embryos into Warm 3 Omni and let the oocytes or embryos remain in the solution for 5 to 10 minutes.

Following warming, oocytes should be held in an incubator at regulated G1+ PLUS, G1+™, G2+™ PLUS or EmbryoGlide™ at 37 °C at 6% CO₂ according to standard laboratory practice.

Specifications

Aseptically filtered
Mouse Embryo Assay, -1°C embryo developed to expanded blastocyst at 96 hours

Bacterial endotoxin (LAL assay) [EU/ml] < 0.5

pH test

Osmolality tested

L-OTC-specific test results are available on the Certificate of Analysis provided with each delivery.

Precautions

The long-term safety of vitrification and/or blastocyst culture on human serum albumin has not been established.

The risk of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife's IVF media, have not been determined.

Discard product if bottle integrity is compromised. Do not use RapidWarm™ Omni if it appears cloudy.

RapidWarm™ Omni contains human serum albumin.

Caution: All blood products should be treated as potentially infectious. Do not use if there is any evidence that this product was derived from a human donor who has had an infection.

Any adverse incident that has occurred in relation to the use of this product should be reported to the manufacturer and the competent authority of the Member State in which the product is marketed.

Not for injection.

Discard the product according to standard clinical practice for medical hazardous waste when the procedure is finished.

Symbols glossary

*Medical devices – Symbols to be used with information to be supplied by the manufacturer (ISO 15223-1,2021).

Symbol Title Reference and symbol ID

Consult instructions for use *Ref no: 5.4.3

Caution *Ref no: 5.4.4

Manufacturer *Ref no: 5.1.1

Date of manufacture *Ref no: 5.1.3

Use-by date *Ref no: 5.1.5

Batch code *Ref no: 5.1.6

Catalogue number *Ref no: 5.1.6

Specimen processing techniques *Ref no: 5.2.2

Single sterile barrier system *Ref no: 5.2.11

Temperature limit *Ref no: 5.3.7

Keep away from sunlight (light)

Contain a medical substance

Medical Device *Ref no: 5.4.7

European conformity (CE) mark identification number as defined in Directive 93/42/EEC 2017/745 Annex V

Annex VI

CE Mark (Conformité Européenne)

CE 2460 CEE/CE (Conformité Européenne)

Caution: Federal (US) law restricts this device to sale by or on the order of a physician or practitioner trained in its use (Rx only).

CHINA

Zh: 玻璃化解冻液RapidWarm™ Omni™²⁺品说明书

产品名称：玻璃化解冻液

适用范围：适用于人类胚胎、卵裂期胚胎和囊胚的玻璃化冷冻后的复苏。

产品型号：型号:RapidWarm™ Omni™, 10124

规格：4.5mL

用途：在培养箱内将温度从-196℃升高至-137℃，同时最高升至-128℃。

RapidWarm™ Omni在包装盒上显示“品名”或“成分”时，该产品分析报告。显示“品名”或“成分”时，该产品分析报告。显示“品名”或“成分”时，该产品分析报告。

注意事项：仅适用于对该设备敏感的患者。

产品特点：清洁白色蛋白。

产品优势：无菌，不含糖。

