SpermGrad™

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en: Indication for use

Medium for gradient sperm separation.

Contraindications

Not known.

Product Description

All SpermGrad™ media are bicarbonate and HEPES buffered media containing silane-coated, colloid silica particles.

SpermGrad is available as:

–SpermGrad REF 10102 (30 mL) and 10099 (125 mL). Media for dilution into desired concentrations. For use with dilution with G-IVF™ PLUS or SpermRinse™.

–SpermGrad REF 10138 (2x30 mL), 10238 (30 mL), 10338 (30 mL), 10139 (2x125 mL), 10239 (125 mL) and 10339 (125 mL) containing: Gradient solutions SpermGrad Upper Layer and SpermGrad Lower Layer. Ready to use solutions for gradient sperm separation.

SpermGrad is antibiotic-free while the media for dilution and sperm wash, such as G-IVF PLUS or SpermRinse, contain gentamicin as an antibacterial agent.

The summary of safety and clinical performance can be found at www.vitrolife.com

SpermGrad is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis provided with each delivery.

Storage instructions and stability

Store dark at +2 to +8°C.

SpermGrad is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis provided with each delivery.

Media bottles should not be stored after opening.

Discard excess media after completion of the procedure.

Directions for use

The product should be used by an IVF professional.

The patient target group is an adult or reproductive-age population that undergoes fertility treatment.

Dilution: SpermGrad REF 10099 and 10102: Mix SpermGrad with G-IVF PLUS or SpermRinse™ in separate tubes to obtain 90% and 45% stock solutions. For 90% stock solution, mix 9.0 mL SpermGrad with 1.0 mL G-IVF PLUS or SpermRinse. For 45% stock solution mix 4.5 mL SpermGrad with 5.5 mL G-IVF PLUS or SpermRinse. Mix the solutions thoroughly. Store in sterile non-toxic tubes or sterile tissue culture flasks. Label with date and refrigerate until use.

SpermGrad REF 10138, 10238, 10338, 10139, 10239 and 10339: Ready to use.

Preparation of sample: Allow the solutions to reach ambient temperature before use.

Layer the gradient in a sterile, conical centrifuge tube marked with patient ID. Pipette 1.5 mL of 90 % solution (SpermGrad Lower Layer) into the tube first and then slowly 1.5 mL of 45 % solution (SpermGrad Upper Layer) on the top of it. Finally, the semen is gently layered on the top. Centrifuge the tube for 20 minutes at 300–600 g remove the two top layers. Transfer the sperm pellet into a sterile conical tube with 5 mL of equilibrated G-IVF PLUS or SpermRinse. Centrifuge for 10 minutes at 300–600 g. Discard the supernatant and repeat the wash. Re-suspend the pellet in 1 mL of equilibrated G-IVF PLUS or SpermRinse. Assess the sample for motility and concentration. Inseminate according to normal laboratory protocol.

SpermGrad REF 10138, 10238, 10338, 10139, 10239 and 10339: Ready to use.

Preparation of sample: Allow the solutions to reach ambient temperature before use.

Layer the gradient in a sterile, conical centrifuge tube marked with patient ID. Pipette 1.5 mL of 90 % solution (SpermGrad Lower Layer) into the tube first and then slowly 1.5 mL of 45 % solution (SpermGrad Upper Layer) on the top of it. Finally, the semen is gently layered on the top. Centrifuge the tube for 20 minutes at 300–600 g remove the two top layers. Transfer the sperm pellet into a sterile conical tube with 5 mL of equilibrated G-IVF PLUS or SpermRinse. Centrifuge for 10 minutes at 300–600 g. Discard the supernatant and repeat the wash. Re-suspend the pellet in 1 mL of equilibrated G-IVF PLUS or SpermRinse. Assess the sample for motility and concentration. Inseminate according to normal laboratory protocol.

SpermGrad REF 10138, 10238, 10338, 10139, 10239 and 10339: Ready to use.

Precautions

Discard product if bottle integrity is compromised. Do not use SpermGrad if it appears cloudy.

Re-use may result in microbiological contamination and/or property changes in the product.

To avoid contamination Vitrolife strongly recommends that media should be opened and used only with aseptic technique.

A very serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Not for injection

Discard the product according to standard clinical practice for medical hazardous waste when the procedure is finished.

Symbols glossary

* Medical devices – Symbols to be used with information to be supplied by the manufacturer (ISO 15223-1:2021).

Symbol	Title	Reference symbol ID
	Consult instructions for use	*Ref no. 5.4.3
	Caution	*Ref no. 5.4.4
	Manufacturer	*Ref no. 5.1.1
	Date of manufacture	*Ref no. 5.1.3
	Use-by date	*Ref no. 5.1.4
	Batch code	*Ref no. 5.1.5
	Catalogue number	*Ref no. 5.1.6
	Sterilized using aseptic processing techniques	*Ref no. 5.2.2
	Single sterile barrier system	*Ref no. 5.2.11
	Do not re-use	*Ref no. 5.4.2
	Temperature limit	*Ref no. 5.3.7
	Keep away from sunlight (light source)	*Ref no. 5.3.2
	Medical Device	*Ref no. 5.7.7

European conformity (CE) mark with Notified Body identification number as defined in Regulation (EU) 2017/745 Annex V

Caution: Federal (US) law restricts this device to sale by or on the order of a physician or practitioner trained in its use (Rx only).

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2024.12 SpermGrad™ © Vitrolife