ThawKit™

Cleave

US

en: Indication for use

Solutions for thawing of frozen procarious oocytes and embryos stored in vitrolife containers.

Contraindications
ThawKit™ Cleave contains gentamicin. Do not use in patients with known hypersensitivity/allergy to the antibiotic.

Product Description

SUPPLEMENTED WITH HSA
ThawKit Cleave contains three solutions for thawing of frozen procarious oocytes and cleavage stage embryos. The solutions consist of MOPS buffered solution containing gentamicin as an antibacterial agent and human serum albumin.

Thawing solution 1 (TS1) and thawing solution 2 (TS2) contain solutions containing decreasing concentrations of sucrose allowing removal of 1,2-propanediol and gradual rehydration of the cells.

Equilibration solution (ES) Buffered solution that does not contain any cryoprotectants used for final rehydration before culture.

Thawing solution 3 (TS3) Buffered solution that does not contain any cryoprotectants used for final rehydration before culture. The quantity of safety and clinical performance can be found at www.vitrolife.com

Storage instructions and stability

Store dark at +2 to +8 °C.

ThawKit Cleave is stable until the expiry date shown on the bottle labels and the LOT specific Certificate of Analysis.

Use by **Two bottles can be used for up to two weeks after first thawing, use aseptic technique and minimize the time outside the refrigerator. Record opening date and time. Discard excess media no later than two weeks after first opening.**

Directions for use

The product shall be used by an IVF professional. Thawing should be performed in an adult or reproductive population that undergoes IVF treatment or fertility preservation.

Thawing should be performed at recommended temperatures and time frames are important. Check patient's identity and prepare all paperwork according to the product.

Label dates, prepare appropriate volumes (0.5-1.0 ml) of TS1, TS2 and ES into respective dishes. Equilibrate embryos in solution 1 for 5 minutes. Always transfer embryos into the TS2 and expose for 5 minutes.

Transfer the embryos to accurately equilibrated culture medium, rins embryos properly and culture according to laboratory procedures.

Specifications

Asepticly filled
Moisture Embryo Assay, 1-cell [% embryos transferred to expanded blastocyst] $\geq$ 80
Bacterial endotoxin (LA, assay) [EU/ml] $<$ 0.5
pH tested
Osmolality tested
LOT-specific test results are available on the Certificate of Analysis provided with each delivery.

Discard product if bottle integrity is compromised. Do not use ThawKit Cleave if it appears cloudy.

ThawKit Cleave contains human serum albumin.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was not found to contain antibodies to HIV, HBs, HCV and HTLV III and non-reactive to HbAg, HCV RNA and HIV RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents and are not.

To avoid contamination Vitrolife strongly recommends that media should be opened and used with aseptic technique.

The risk of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife V media, have not been determined and are uncertain.

Any serious incident that has occurred in relation to the product should be reported to the manufacturer and the user and/or patient is established.

Not for injection.
Discard product according to standard clinical practice for medical hazardous waste when the procedure is finished.

Symbol Legend
ThawKit Cleave - Symbols to be used with information to be supplied by the manufacturer (ISO 15222-1:2001).

Symbol	Title	Referenced symbol ID
	Consult instructions for use	Ref no: 5.4.3
	Caution	Ref no: 5.4.4
	Manufacturer	Ref no: 5.4.1
	Date of manufacture	Ref no: 5.1.3
	Use by date	Ref no: 5.1.4
	Batch code	Ref no: 5.1.5
	Catalogue number	Ref no: 5.1.6
	Sterilized using aseptic technique	Ref no: 5.2.2
	Single use system	Ref no: 5.2.1
	Temperature	Ref no: 5.3.1
	Keep away from sunlight (light protection)	Ref no: 5.3.2
	Contains a medicinal substance	Ref no: 5.4.7
	Contains human blood or plasma	Ref no: 5.4.6
	Medical Device	Ref no: 5.7.7

Symbol	Title	Referenced symbol ID
	CE Mark (Conformity European)	Ref no: 5.7.7

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Vitrolife Sweden AB, Gustaf Werners gata 2, SE-421 32 Västra Frölunda, Sweden
Vitrolife Inc., 3601 South Inca Street, Englewood CO 80110, USA
Vitrolife Pty. Ltd., Level 10, 68 Pitt Street, Sydney NSW 2000, Australia
Manufactured by Vitrolife Sweden AB