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Active	Name	Type	Creator	Date
<input checked="" type="checkbox"/>	Model name	Additive	ADMIN	2014-10-28

Joonis 1: Nupp Import (Impordi) lehel Settings (Seaded) vahelehel Models (Mudelid).

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The screenshot displays the 'Models' tab in the KIDScore software. On the left, there is a table listing models:

Active	Name	Type	Creator	Date
<input checked="" type="checkbox"/>	KIDScoreD3 v1.4	Imported	Vitrolife	2024-09-13
<input checked="" type="checkbox"/>	KIDScoreD5 v3.3	Imported	Vitrolife	2024-09-13

Below this table are 'Import' and 'Export' buttons. A 'Selected model' section shows a list of variables with 'Min' and 'Max' values:

Variable	Min	Max
Information - NOT3PN	-	-
Information - 9Pnf	-	-
Information - t2	-	-
Information - t3	-	-
Information - t4	-	-
Information - t5	-	-
Information - t8	-	-
Information - Cells 66h	-	-

The main area shows the details for 'KIDScore D3 v1.4'. The 'Model Name' is 'KIDScoreD3 v1.4', 'Model Type' is 'Imported', and 'Creator' is 'Vitrolife'. The 'Model Description' states: 'KIDScore D3 is defined by Vitrolife A/S based on the knowledge and experience extracted from our available KID data (please see the use manual for the EmbryoViewer software for a definition of KID data). The model focuses on which embryos to avoid rather than which embryos to select. It is thus a model which is based on avoidance criteria rather than selection criteria. The model will apply the avoidance criteria to the embryo and...'. Below this is a 'Model Definition' table:

Variable	Description	Min	Max	Classification
NOT3PN	Info			
9Pnf	Info			
t2	Info			
t3	Info			
t4	Info			
t5	Info			
t8	Info			
Cells 66h	Info			

At the bottom right, there is a license agreement section for 'Vitrolife A/S End-User License Agreement for KIDScore D3 Model'. It includes the Vitrolife logo and contact information: 'KIDScore D3 Software version 1 REF 16531 VERSION 1.4.0.29558'. The license agreement text states: 'The installation and use of this KIDScore D3 model (the "Model") shall be subject to the terms and conditions stated below. By clicking the "I agree" button or by installing or otherwise using the Model you have conclusively accepted to be bound by all of these terms and conditions. If you do not agree to all of the terms and conditions please do not install or use the Model. All rights in the Model belong to Vitrolife A/S ("Vitrolife"). Vitrolife grants you a non-exclusive, non-transferable, and non-sub-licensable license to install the Model and use it solely together with your EmbryoViewer and ES Server software with a valid license from Vitrolife. You are not granted any other rights or license with respect to the Model. Without limiting the foregoing you shall not copy, modify, decompile, reverse engineer, disassemble, or convert the Model or assign, transfer, sell, rent or lease the Model to any third party. Any actions, use, copying or distribution of the Model not authorized under these terms of use shall automatically terminate your rights hereunder. VITROLIFE DISCLAIMS ALL WARRANTIES EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTIES OF FITNESS FOR A PARTICULAR PURPOSE, MERCHANTABILITY OR NONINFRINGEMENT OF THIRD-PARTY RIGHTS. VITROLIFE SHALL NOT UNDER ANY CIRCUMSTANCES BE LIABLE TO YOU OR ANY THIRD PARTY FOR ANY LOSS OF PROFITS, LOSS OF USE, INTERRUPTION OF BUSINESS, OR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY KIND. This end-user license agreement shall be governed and interpreted solely in...'

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
The screenshot displays the EmbryoViewer interface. On the left, a large circular time-lapse image shows an embryo at 111.9h. Below the image is a timeline with colored segments (green, yellow, blue) and a red signal graph. On the right, a table lists developmental stages with their corresponding values and progress percentages.

Anno. Tool	Value	Well	Dec.	Progress
PN	2	AA-1		100 %
t2	28.5 h	AA-2		0 %
t3	30.9 h	AA-3		0 %
t4	31.2 h	AA-4		0 %
t5	42.8 h	AA-5		0 %
tB	77.6 h	AA-6		0 %
ICM	A	AA-7		0 %
TE	B	AA-8		0 %
		AA-9		0 %
		AA-10		0 %
		AA-11		0 %
		AA-12		0 %
		AA-13		0 %
		AA-14		0 %
		AA-15		0 %
		AA-16		0 %

Below the table, there is a 'Comments' field, 'Prev' and 'Next' navigation buttons, a 'Strategy' dropdown menu set to 'KIDScore D5 v3 Annotation', and radio buttons for 'Step by' with 'Variable' selected and 'Well' unselected.

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Well	Dec.	Current score	NOT2PN	t2	t3	t4	t5	tB	ICM	TE	Last stage	Morph. grade	Last image	Saved score
AB-1		6.1	●	25.5	38.3	39.6	50.5	109.7	B	B		B		
AB-2		6.6	●	25.8	38.3	40.2	53.6	103.7	B	B		EB		
AB-3		8.4	●	27.7	40.5	41.3	53.9	103.9	B	A		EB		
AB-4		6.2	●	25.8	39.1	39.2	55.5	106.6	B	B		EB		

Current Model
 KIDScoreD5 v3
 Created 2018-11-01 by Vitrolife

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Well	Dec.	Current score	UNEVEN2	UNEVEN4	MN2	MN4	Frag-2	Coll. Count	Last stage	Morph. grade	Last image	Saved score
AB-1		NA	●	●	●	●	5.0	0.0		B		6.1
AB-2		NA	●	●	●	●	5.0	1.0		EB		6.6
AB-3		NA	●	●	●	●	20.0	0.0		EB		8.4
AB-4		NA	●	●	●	●	10.0	2.0		EB		6.2

Current Model
 Info example
 Created 2019-06-04 by ADMIN

Saved Model
 KIDScoreD5 v3
 Saved 2019-07-02 15:57:22 by ADMIN

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©2024 Vitrolife A/S. Kõik õigused kaitstud.
VITROLIFE A/S, TAANI