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Active	Name	Type	Creator	Date
<input checked="" type="checkbox"/>	Model name	Additive	ADMIN	2014-10-28

Slika 1: Dugme „Import“ (Uvezi) na stranici „Settings“ (Podešavanja), na kartici „Models“ (Modeli).

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The screenshot displays the 'Models' tab in the KIDScore software. It features a table of models, a 'Selected model' section, and a detailed view of the 'KIDScore D3 v1.3' model. The detailed view includes the model name, type, creator, and description, along with a table of variables and their classifications.

Active	Name	Type	Creator	Date
<input checked="" type="checkbox"/>	KIDScoreD5 v3.2	Imported	Vitrolife	2019-10-01
<input checked="" type="checkbox"/>	KIDScoreD3 v1.3	Imported	Vitrolife	2022-03-16

Selected model

Variable	Min	Max
Information - NOT2PN	-	-
Information - IPNF	-	-
Information - t2	-	-
Information - t3	-	-
Information - t4	-	-
Information - t5	-	-
Information - t8	-	-
Information - cells66	-	-

Model Definition

Variable	Description	Min	Max	Classification
NOT2PN	Info			
IPNF	Info			
t2	Info			
t3	Info			
t4	Info			
t5	Info			
t8	Info			
cells66	Info			

Model Description

KIDScore D3 is defined by Vitrolife A/S based on the knowledge and experience extracted from our available KID data (please see the use manual for the EmbryoViewer software for a definition of KID data).

The model focuses on which embryos to avoid rather than which embryos to select. It is thus a model which is based on avoidance criteria rather than selection criteria. The model will apply the avoidance criteria to the embryos and...

Model Provided By:

Vitrolife

Vitrolife A/S End-User License Agreement for KIDScore D3 Model.

The installation and use of this KIDScore D3 model (the "Model") shall be subject to the terms and conditions stated below. By clicking the "I agree" button or by installing or otherwise using the Model you have conclusively accepted to be bound by all of these terms and conditions. If you do not agree to all of the terms and conditions please do not install or use the Model.

All rights in the Model belong to Vitrolife A/S ("Vitrolife"). Vitrolife grants you a non-exclusive, non-transferable, and non-sub-licensable license to install the Model and use it solely together with your EmbryoViewer and ES Server software with a valid license from Vitrolife.

You are not granted any other rights or license with respect to the Model. Without limiting the foregoing you shall not copy, modify, decompile, reverse engineer, disassemble, or convert the Model or assign, transfer, sell, rent or lease the Model to any third party. Any actions, use, copying or distribution of the Model not authorized under these terms of use shall automatically terminate your rights hereunder.

VITROLIFE DISCLAIMS ALL WARRANTIES EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTIES OF FITNESS FOR A PARTICULAR PURPOSE, MERCHANTABILITY OR NONINFRINGEMENT OF THIRD-PARTY RIGHTS. VITROLIFE SHALL NOT UNDER ANY CIRCUMSTANCES BE LIABLE TO YOU OR ANY THIRD PARTY FOR ANY LOSS OF PROFITS, LOSS OF USE, INTERRUPTION OF BUSINESS, OR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES OF ANY KIND.

This end-user license agreement shall be governed and interpreted solely in...

KIDScore D3
Software version 1
REF: 15531
VERSION 1.3.0.27527

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
The screenshot displays the KIDScore software interface. On the left, there is a large circular image of an embryo. Below it, a timeline shows various stages with colored bars (green, yellow, blue) and a red signal graph. The top right corner shows 'Embryo ID: AA1'. The main right panel contains a table with columns: Anno. Tool, Value, Well, Dec., and Progress.

Anno. Tool	Value	Well	Dec.	Progress
PN	2	AA-1		100 %
t2	28.5 h	AA-2		0 %
t3	30.9 h	AA-3		0 %
t4	31.2 h	AA-4		0 %
t5	42.8 h	AA-5		0 %
tB	77.6 h	AA-6		0 %
ICM	A	AA-7		0 %
TE	B	AA-8		0 %
		AA-9		0 %
		AA-10		0 %
		AA-11		0 %
		AA-12		0 %
		AA-13		0 %
		AA-14		0 %
		AA-15		0 %
		AA-16		0 %

Below the table, there is a 'Comments' field, 'Prev' and 'Next' buttons, a 'Strategy' dropdown menu set to 'KIDScore D5 v3 Annotation', and radio buttons for 'Step by' with 'Variable' selected and 'Well' unselected.

3. Na stranici **View Running** (Prikaz aktivnog) ili **View All Slides** (Prikaz svih slajdova) izaberite označenu EmbryoSlide posudu za kulturu. Zatim kliknite na **Compare & Select** (Upoređivanje i odabir) na navigacionom panelu.

4. Sa padajuće liste na stranici **Compare & Select** (Upoređivanje i odabir) izaberite željeni model i primenite ga na embrione u aktuелnoj EmbryoSlide posudi za kulturu.

Model sada dodeljuje numeričku ocenu svakom embrionu. Embrioni sa najvišim ocenama imaju statistički najveću šansu za implantaciju. Međutim, imajte u vidu da parametri koji nisu uključeni u model takođe mogu da utiču na kvalitet embriona.

Well	Dec.	Current score	NOT2PN	t2	t3	t4	t5	tB	ICM	TE	Last stage	Morph. grade	Last image	Saved score
AB-1		6.1	●	25.5	38.3	39.6	50.5	109.7	B	B		B		
AB-2		6.6	●	25.8	38.3	40.2	53.6	103.7	B	B		EB		
AB-3		8.4	●	27.7	40.5	41.3	53.9	103.9	B	A		EB		
AB-4		6.2	●	25.8	39.1	39.2	55.5	106.6	B	B		EB		

Current Model
 KIDScoreD5 v3
 Created 2018-11-01 by Vitrolife

Saved Model
 No saved model

Transfer Info
 Transfer Date: 2019-07-02

View All Patient Embryos

5. Kliknite na **Save** (Sačuvaj).
6. Možete da primenite dodatni model na embrione kako biste uzeli u obzir i druge parametre koji nisu obuhvaćeni KIDScore modelom. Da biste to uradili:
- Izaberite korisnički definisan model i ocenite embrione korišćenjem tog modela.
 - Proverite i uporedite rezultate. Ocene koje je dodelio KIDScore model u prethodnim koracima pojavljuju se u koloni pod nazivom **Saved score** (Sačuvana ocena).

Primeri korisnički definisanih parametara

Well	Dec.	Current score	UNEVEN2	UNEVEN4	MN2	MN4	Frag-2	Coll. Count	Last stage	Morph. grade	Last image	Saved score
AB-1		NA	●	●	●	●	5.0	0.0		B		6.1
AB-2		NA	●	●	●	●	5.0	1.0		EB		6.6
AB-3		NA	●	●	●	●	20.0	0.0		EB		8.4
AB-4		NA	●	●	●	●	10.0	2.0		EB		6.2

Min
Max
Weight

Current Model
 Info example
 Created 2019-06-04 by ADMIN

Saved Model
 KIDScoreD5 v3
 Saved 2019-07-02 15:57:22 by ADMIN

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©2022 Vitrolife A/S. Sva prava zadržana.
VITROLIFE A/S, DANSKA