

# EMBRYO CULTURE GX-TL

## Day 0



Gx-TL™



OVOIL™

Prepare micro-droplet culture dishes with 25 µL droplets of Gx-TL for washing and for culture. Cover with OVOIL and pre-equilibrate at

**37°C 6 % CO<sub>2</sub>  
overnight**



Gx-MOPS™ PLUS

Warm Gx-MOPS PLUS (for fertilisation assessment) in rinsed tightly capped tubes in a warming incubator **without CO<sub>2</sub>\*** at

**37°C overnight**

\* An adequately calibrated warming block can be used for tubes instead of a warming incubator.


**Ensure that the denudation and washing procedures are performed at 37°C**

## Day 1

### 1. Fertilisation assessment

For inseminated oocytes, transfer the oocytes to a centre well dish with pre-warmed Gx-MOPS PLUS. If denudation and fertilisation assessment can be performed within 2 minutes, Gx-IVF can be used instead of Gx-MOPS PLUS. Remove cumulus and corona cells from oocytes using a denudation pipette and assess fertilisation at

**37°C**



For ICSI oocytes, assess fertilisation in the Gx-TL micro-droplet culture dish.


### 2. Culture

Wash the zygotes extensively in the Gx-TL micro-droplet dishes prepared on Day 0 and transfer the zygotes to 25 µL Gx-TL culture droplets covered with OVOIL. Culture to the blastocyst stage at



**37°C 6 % CO<sub>2</sub>**

If your clinic is located at a higher altitude than sea level, CO<sub>2</sub> percentage should be increased, see graph below.




**We recommend Gx-MOPS PLUS for washing and assessment of embryos before transfer to Gx-TL or EmbryoGlue®**

## Day 5

### In the morning of day 5

Assess embryo cleavage and move the blastocysts selected for transfer to the equilibrated transfer dish and leave at



**37°C 6 % CO<sub>2</sub>  
until transfer**

**For blastocyst transfer, see separate Embryo transfer protocol**