

DISH LOADING

Primo Vision time-lapse system

Dish preparation


9- or 16 well dish

Use a flexible pipette with a diameter of 100 μm .

Place the tip of the pipette at the bottom of the wells when filling with media.



Fill each well with pre-warmed media to the top of the well but not exceeding the well.



Repeat this step one by one for each well of the dish.



When each well is filled use the remaining medium in the tip to connect the wells.

Then enlarge the drop to
80 μl : 9-well dish
150 μl : 16 well dish

Spread the drop to reach a diameter of approximately

9 mm: 9-well dish
13 mm: 16-well dish



Cover with 3-4 ml of oil.
Use a stereo microscope to ensure sufficient overlay.


Equilibrate for 16 hours

Removing bubbles

After equilibration, remove occasional bubbles by dropping the dish gently to a hard surface. Remove the lid before dropping the dish to avoid oil spilling to the lid.

This will bring occasional bubbles out of the wells and into the media droplet.


Remove bubbles from the media droplet using a standard pipette.



Loading embryos

Load embryos into each well using a micropipette.

Release embryos at the top of the well and allow them to settle at the bottom of the wells.



Place the dish under the EVO+ microscope in your incubator.

Place the dish under the EVO+ microscope inside the incubator