Embryo culture Gx-TL

Day 0



Prepare micro-droplet culture dishes with 25 µL droplets of Gx-TL for washing and for culture. Cover with OVOIL and pre-equilibrate at

> 37°C 6 % CO₂ overnight



Warm Gx-MOPS PLUS (for fertilisation assessment) in rinsed tightly capped tubes in a warming incubator without CO₂* at

37°C overnight

*An adequately calibrated warming block can be used for tubes instead of a warming incubator.

Ensure that the denudation and washing procedures are performed at 37°C

Day 1

1. Fertilisation assessment

For inseminated oocytes, transfer the oocytes to a centre well dish with pre-warmed Gx-MOPS PLUS. If denudation and fertilisation assessment can be performed within 2 minutes, Gx-IVF can be used instead of Gx-MOPS PLUS. Remove cumulus and corona cells from oocytes using a denudation pipette and assess fertilisation at



For ICSI oocytes, assess fertilisation in the Gx-TL micro-droplet culture dish.

2. Culture


Wash the zygotes extensively in the Gx-TL micro-droplet dishes prepared on Day 0 and transfer the zygotes to 25 µL Gx-TL culture droplets covered with OVOIL. Culture to the blastocyst stage at

37°C 6 % CO₂

If your clinic is located at a higher altitude than sea level, CO_2 percentage should be increased, see graph below.

12.0 %	+					
11.0 %	-					/
10.0 %	-					
9.0%	-					
8.0 %	-					
7.0 %	-					
6.0%		-				
5.0 %						T
() m	1000 m	2000 m	3000 m	4000 m	5000 m

We recommend Gx-MOPS PLUS for washing and assessment of embryos before transfer to Gx-TL or EmbryoGlue®



Vitrolife 🔨