

pH validation of culture medium in EmbryoScope Flex


For optimal embryo development, it is necessary to maintain a proper pH in the growth medium. The optimal pH may differ for different media, but a pH between 7.2 and 7.4 is frequently recommended. For long term incubation of human cleavage stage embryos most IVF laboratories use a culture medium that is based on a bicarbonate buffer system and requires a controlled atmosphere with elevated carbon dioxide (5 – 8%).

The control of the proper pH is mandatory in quality assurance. There are two possibilities to achieve this goal: 1) measuring pH and 2) measuring CO₂. Measuring pH is the only way to detect slight changes between different lots of the same culture media at exactly the same CO₂ concentration.

Here we present a combined method to measure pH and CO₂ in the EmbryoScope Flex time-lapse system. Once the correlation between pH and CO₂ of a given media lot has been established, further validations can be done by measuring only CO₂ until a new lot of media is used, requiring a new pH/CO₂ correlation measurement.


Carbon dioxide and pH of IVF culture medium

The pH of the medium is controlled by the carbon dioxide concentration in the incubation chamber. The figure below shows the theoretical relationship between CO₂ concentration in an incubator (in % atm) and the resulting equilibrium pH as a continuous blue curve. The green points are experimental values from *D. Gardner and M. Lane (2000) Embryo culture systems in Handbook of In Vitro Fertilization, 2nd Ed, CRC Press p. 232-234.*



Calculation based on a medium with a total alkalinity of 27.4 mmol/kg at 37°C, CO₂SYS, van Heuven et al 2009: http://cdiac.ornl.gov/ftp/co2sys/CO2SYS_calc_MATLAB/

Carbon dioxide in the incubator is in equilibrium with dissolved carbon dioxide, carbonic acid, bicarbonate and carbonate:



The pH of the media at physiological pH can be approximated by the Henderson-Hasselbach equation:

$$\text{pH} = \text{pK}_a + \log_{10}([\text{HCO}_3^-]/[\text{CO}_2])$$

The relationship between carbon dioxide concentration within an incubator and the resulting pH in a given medium shows a remarkable resilience to pH change even after substantial change in CO₂ concentration. Large changes in CO₂ thus only cause minor changes in pH. We recommend a stable CO₂ measuring device that gives accurate results in the desired range of 5-8%.

Validation of medium pH for a specific carbon dioxide concentration

It is important to validate the pH of the media being used on a regular basis, to make sure the pH follows recommendations for the specific media. This validation can be performed by measuring the pH of a media sample in a culture dish placed inside the incubation chamber for at least 24 hrs. For this purpose, we recommend to fill the EmbryoSlide Flex culture dish with 210µL into the reservoir and the two flushing wells with 35 µL of culture medium under a 1.2mL culture oil overlay. This will facilitate a sample size of at least 150µl for analysis. It is very important to add culture oil in order to avoid evaporation as the environment inside the EmbryoScope Flex time-lapse incubator is NOT humidified.

The method is described on the next page and requires a blood gas analyser.

How to measure pH of culture medium in the EmbryoScope Flex incubator.

1. In the morning before opening the EmbryoScope Flex incubator measure the CO₂ concentration in the incubation chamber using a calibrated CO₂ measuring device as described in the User Manual. Calibrate if necessary.
2. Fill the EmbryoSlide Flex culture dish with 210µL of culture media into the reservoir and eventual 35µL into each two flushing wells.
3. Overlay with 1.2mL of culture oil.
4. Place the dish in the EmbryoScope Flex incubator and let equilibrate for 24 hours.
5. After 24 hours remove the EmbryoSlide Flex culture dish from the incubation chamber and according to the specification of the blood gas analyser immediately pipette the appropriate amount of medium from the reservoir and flushing wells. Avoid aspirating oil.
6. Transfer the sample to a calibrated blood gas analyser and take the pH measurement immediately.
7. Until a new lot of culture medium is used, regular CO₂ measurements are sufficient as the correlating pH value is known.
8. If the pH measurement is lower than described in your laboratory standards, it means, the CO₂ concentration is too high. To compensate for this, change the set-point on the incubator to a lower level and repeat the pH measurement procedure as just described with a new dish with fresh medium and equilibrate 24 hours before a new pH measurement is performed.

Note

When measuring and interpreting media pH it is important to remember that:

1. The pH-CO₂ relationship is altitude dependant. That means that at altitudes higher than sea level, CO₂ % must be higher to obtain the same pH.
2. Molecular diffusion through the oil and media layer to reach a stable equilibrium takes several hours. It is thus essential that the media sample is *equilibrated for at least 24 hrs* before a stable reliable media pH can be measured (Einstein-Smoluchowski equation: $t=s^2/2D$).
3. When handling small media samples care must be taken to avoid temperature changes (see figure to the right) and CO₂ diffusion into plastware. Pre-equilibration of plastware under CO₂ atmosphere may be necessary. The figure (right) shows the moderate effect of reducing media temperature to room temperature while maintaining 5% CO₂. Cooling the media to RT will decrease pH by 0.07 which must be corrected for.

