

Software Release Note

kMap Software Release Notes: 23Jun2025

UI Release Version	Algorithm Version	Update
0.1.5	0.2.24 (unchanged)	Improved session handling to prevent unexpected timeouts.
		Automatic removal of Bead Arrays from the system once all associated samples have been deleted from a case.
		Complete display of all Samples and Files within the user interface for each case.
		Resolved intermittent loading problems preventing Bead Arrays data from displaying in some cases.
		Upgraded database engine for improved security and performance.
		Enhanced kMap Uploader tool will be made available via link in the software banner with installation instructions: Improved error handling, stricter authentication, better text encoding support, and masked passwords in logs.