

Software Release Note

kMap Software Release Notes: 17Mar2025

UI Release Version	Algorithm Version	Update
0.1.4	0.2.24 (unchanged)	Added security feature to enable Brute Force Protection for Group accounts.
		Added security feature to enable Multi-Factor Authentication for Group accounts.
		Expanded gene list, now covering over 40,000 HGNC gene names used as part of the auto-complete selector. <i>Note: Gene list is provided to support users, and users are advised to confirm their regions of interest.</i>
		Improved usability with text sizing and sample separation.
		Improved responsiveness of sample processing pipeline.
		Improved consistency of the Trio icon status shown between the analysis and the generated report when analysing multiple target regions in a case.
		Software-wide usability and logging improvements.